

# Current issues in *Aspergillus* quantitative real time PCR (RTqPCR) standardization



**Monica Herrera MD and Brian Wickes PhD**

*Dept of Microbiology and Immunology, The University of Texas Health Science Center at San Antonio, TX*

**IAAM**



INVASIVE ASPERGILLOSIS ANIMAL MODELS

S ♦ A ♦ C ♦ M ♦ M

San Antonio Center for Medical Mycology

# Background

---

- Ubiquitous nature makes exposure difficult to avoid.
- IA life- threatening nature makes accurate diagnosis and early detection crucial.
- Current methods need improved sensitivity
- RTqPCR offers a rapid and sensitive method that needs to be calibrated in order to be tested between institutions.
- Lack of validated RTqPCR data remains a limiting factor for standardization and calibration between laboratories

# RTqPCR Assay Variables

## Instrumentation

- Various instruments different thresholds settings.
- Need to standardize settings between laboratories.

## Reagents

- Taq polymerases, buffers probes, primers

## Target

- multi copy rDNA
- FKS1 (1 copy)
- Conservation (primer design)
- amplification efficiency

## Template preparation

- specimen type
- tissue vs. fluid
- DNA isolation methodology
  - yield
  - purity

## Cycling parameters

- Ct values for positive samples should be validated by a standard curve. Limiting threshold should be established.


## Data interpretation

- Allowable limit of detection Ct #of replicates to detect interference effect with 95% confidence power

# RTqPCR Targets: Multicopy


## rDNA Organization of the *Aspergillus fumigatus* Genome

29.4 Megabases  
9,926 ORFs




→ = 1 repeat

AF293 has 38 repeats




# Significance of rDNA in Fungal Detection and Identification




- 1) rDNA subunits are highly conserved
  - ✓ *ITS, and D1/D2 regions are species-specific – variable region*
  - ✓ *Allows universal PCR primer or probe sites in 18s and 28s regions*
- 2) Multicopy nature enhances PCR sensitivity 38-100x

# RTqPCR Targets: Single copy vs. multiple copies

Af293  
*FKS1* vs. *ARG4* and *pyrG*



Af293  
18s rDNA vs *FKS1*





## *A. fumigatus* multicopy rDNA Genes Vary in Copy Number

| STRAIN  | 18s rDNA copy # |
|---------|-----------------|
| Af293   | 38 +/- 0.01     |
| WSA-450 | 42 +/- 0.07     |
| WSA-172 | 46 +/- 0.03     |
| WSA-446 | 47 +/- 0.01     |
| WSA-445 | 49 +/- 0.06     |
| WSA-271 | 49 +/- 0.05     |
| WSA-270 | 53 +/- 0.01     |
| WSA-621 | 70 +/- 0.03     |
| WSA-419 | 91 +/- 0.03     |

Range: 38-91, Avg: 54

Af293 genome=35, QRT-PCR=38

# Aspergillus DNA isolation/ quantitation for standards construction






## A. *Fumigatus* DNA quantification gel vs. nanodrop

### DNA Concentration by gel quantification

| Sample # | gel ng/ul        | nanodrop ng/ul  |
|----------|------------------|-----------------|
| AF293#1  | 9.2              | 447             |
| AF293#2  | 9.8              | 656             |
| AF293#3  | 10.8             | 422             |
| AF293#4  | 8.2              | 364             |
| AF293#5  | 10               | 1723            |
| AF293#6  | 12               | 2861            |
| AF293#7  | 11.2             | 1496            |
|          | <b>AVG 50.86</b> | <b>AVG 1138</b> |

# Primer Efficiencies for FKS1 and 18srDNA genes




## Primer Efficiencies


|                |                |
|----------------|----------------|
| <b>FKS1</b>    | <b>18srDNA</b> |
| <b>0.87955</b> | <b>0.85172</b> |

Primer efficiencies are important because they are a measure of the sensitivity and efficiency of our assay.

# Aspergillus DNA isolation from Biological fluids, tissue





# Sample processing for fungal tissue DNA extraction





RTqPCR Results of GP infected aspergillus lung homogenate during bead beating at different time intervals

# Sample processing for fungal blood DNA extraction



RTqPCR results of GP infected aspergillus blood samples processed through a **column vs. automated** protocol

# Sample analysis in RTqPCR of fungal samples



A standard curve is fundamental to establish valid thresholds for calling a sample a true **positive** or **negative**.

# Summary and Conclusions

- Calibration of RT-qPCR using standard curves enhances accuracy.
  - ✓ Standard curves can identify technical errors
  - ✓ Standard curves are essential for interlaboratory agreement
- Template extraction methods greatly affect results.
  - ✓ DNA standards should be high purity
  - ✓ Automated extraction may enhance sensitivity
- Instrumentation parameters also greatly affect results.
  - ✓ Ct values can be arbitrarily adjusted

# Acknowledgments

Wickes Lab, Dept of Microbiology, UTHSCSA

Jianmin Fu, PhD

Patterson Lab, Dept of Medicine, UTHSCSA

Tom Patterson, MD

Ana Vallor, PhD

Laura Najvar

Rick Kirkpatrick

Dr Caliendo MD, PhD, Emory University

Jessica Ingersoll

AsTec

IAAM